

Interview OphtalmoTV : Ophtalmic PerfeXion et OphtalBoX

Dr Patrice VOTAN interroge Dr Jean-Philippe COLLIOT sur les nouvelles lentilles pour presbytes astigmatas Ophtalmic PerfeXion HR RX Toric Prog et le meuble connecté OphtalBoX permettant de stocker jusqu'à 1000 lentilles d'essai Ophtalmic.

http://www.ophtalmo.tv/La-lentille-Ophtalmic-PerfeXion-et-l-OPHTALBOX-Interview-du-Dr-Jean-Philippe-COLLIOT-par-le-Dr-Patrice-VOTAN-OphtalmoTV_v2134.html



The screenshot displays the OphtalmoTV website interface. At the top left is the OphtalmoTV logo, and at the top right is a search bar labeled "Recherche". Below the logo is a navigation menu with the following items: Accueil, Catégories, BANDES-ANNONCES, Congrès, Symposiums, Qui sommes-nous ?, and Nous écrire. The main content area features a large title: "La lentille Ophtalmic PerfeXion et l'OPHTALBOX : Interview du Dr Jean-Philippe COLLIOT par le Dr Patrice VOTAN - OphtalmoTV". To the right of the title are two small icons: a double-headed arrow and a share icon. Below the title is a metadata line: "Par : Dr Jean-Philippe COLLIOT - Chantilly & Dr Patrice VOTAN - OphtalmoTV | Durée : 8min 40sec | Postée : 02/09/2016 | Chaîne : Contactologie". The video player itself shows a man with a beard (Dr Jean-Philippe COLLIOT) speaking, with a small "Ophtalmic" logo in the bottom left corner of the video frame. To the right of the video player are social media icons for Twitter and Facebook. The background of the video player is a light blue color with a grid pattern.